






Gynecologic Testing Options

Collection Kit	Sample Type	Pap	HPV	HPV 16 18/45 Genotype Assay	HSV 1 & 2*	CT/NG Assay*	Vaginitis Panel* (Note changes below)	<i>Mycoplasma genitalium</i> *
ThinPrep® (PreservCyt®)	 Vaginal or Cervical/endocervical	✓	✓	✓	✓	✓		COMING SOON
Aptima® Multitest Swab Collection Kit	 Clinician or patient collected vaginal swab, Clinician collected male urethral swab					✓	✓	
Aptima® Urine Specimen Collection Kit	 Urine (male or female)					✓		

* HSV/STI testing must be ordered in conjunction with other testing on a single specimen, as it cannot be performed on material previously processed.



Vaginitis Panel collection device has changed!! No more ampoules to break. Same Aptima® orange collection swab as for CT/NG testing.

Our Vaginitis Panel now reports each of the following:

Bacterial vaginosis: Targets *Gardnerella vaginalis*, *Atopobium vaginae* & *Lactobacillus species* and returns a simple qualitative BV positive or negative result.

CV/TV: Differentiates between *Candida species* group (composed of *C. albicans*, *C. tropicalis*, *C. parapsilosis*, *C. dubliniensis*) and *C. glabrata* to aid in determining the most appropriate antifungal therapy for patients, and confirms the presence of a *Trichomonas vaginalis* infection (recommended by the CDC in all women seeking treatment for vaginal discharge). All are reported with either a positive or negative result.

- Sexually Transmitted Diseases Treatment Guidelines, 2015. <http://www.cdc.gov/std/tg2015/trichomoniasis>. Reviewed June 4, 2015. Accessed August 19, 2019).

Vista Pathology Laboratory, LLC • 1032 E. Jackson Street • Medford, OR 97504 • Phone (541) 770-4559 • Fax (541) 770-4511